

GelA-GTA Hydrogel Formation

15% GelA wt./vol

2.5 mL 15%
GelA Soln

Various GTA
Solns.


GelA-GTA
mol/mol

Cross-Linked:
60 min- 23°C

Volume

1.00 mL 0.50 mL 0.25 mL

Heated to 65 C°



Stirred



Concentration

2.50% GTA



1.50% GTA

0.50% GTA

	1	2	3
A	15	30	60
B	25	50	100
C	75	150	300

GTA Solution Conditions

Cross-Linking Density



Molar Ratio GelA-GTA

